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**EN****ATTENTION!**

- When purchasing the device check the following:
  - the complete set of delivery (page 31);
  - compliance of serial numbers in this Operating Manual (page 91) with the numbers on the device;
  - availability of sale date and stamp of the seller in the Warranty Card (page 91).
- **Do not switch on the device if it was kept at temperatures below 0 °C!** The device should be kept in conditions of room temperature for at least 2 hours.
- Before operating this device, read this Manual carefully and keep it during the entire useful life of the device.

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Dear customer!

We congratulate you on purchasing the RIKTA® – ESMIL® (2A) device and we express confidence that it will help to maintain your health for many years!

**INVISIBLE LASER RADIATION!**

**DO NOT STARE INTO THE BEAM AND DO NOT MAKE  
DIRECT OBSERVATIONS USING OPTICAL INSTRUMENTS!**

**CLASS 1M LASER EQUIPMENT**

## TABLE OF CONTENTS

<b>RU</b> .....	3
<b>EN</b>	
INTENDED USE .....	26
GENERAL INFORMATION .....	27
SAFETY PRECAUTIONS .....	28
DELIVERY SET .....	31
DEVICE .....	32
PREPARATION FOR OPERATION .....	34
DEVICE ACTIVATION .....	34
DEVICE SERVICEABILITY CHECK (TESTING) .....	35
CONDUCTING THE THERAPY SESSION .....	36
TENS-THERAPY .....	38
BUILT-IN BATTERY MONITORING .....	40
BATTERY RECHARGING .....	40
DEVICE TURNING OFF .....	41
TECHNICAL SPECIFICATIONS .....	42
MAINTENANCE .....	43
POSSIBLE MALFUNCTIONS AND TROUBLESHOOTING .....	44
PACKAGE, STORAGE AND TRANSPORTATION .....	45
MANUFACTURER'S WARRANTY .....	46
<b>FR</b> .....	47
<b>DE</b> .....	69
ACCEPTANCE CERTIFICATE .....	91
GUARANTEE CARD .....	91

### DESIGNATIONS WHICH ARE USED IN OPERATING MANUAL

warning about the situation,  
in which it is possible  
to cause harm to health

mandatory action symbol

notes and tips

RIKTA® – ESMIL® (2A) device is manufactured at the enterprise JSC «MILTA – PCP GIT», whose quality management system complies with the requirements of international standard ISO 13485:2016.

The device is entered in the State register of medical devices of the Russian Federation; its registration certificate No. P3H 2013/1003 has no expiration date.

The device is certified in the Russian Federation.

The Certificate of conformity No. РОСС RU.ИМ02.Н17822 is valid until 26.08.2019.

## INTENDED USE

The portable electro-photo-magnetic-infrared laser therapy RIKTA® – ESMIL® (2A) device is intended for pain relief and non-pharmacological, painless and noninvasive (without breaking the skin integument) treatment of wide range of human diseases. The device can be used simultaneously with medicinal methods of treatment and other methods of therapy.

The device is intended for use in polyclinics, clinics, hospitals, treatment rooms, as well as at home after consultation with a laser therapy specialist.

RIKTA® – ESMIL® (2A) device provides simultaneous penetrative action on the organism of low-intensity laser radiation of the near-infrared range in combination with a constant magnetic field, visible red light and broadband infrared incoherent radiation (**photo-magnetic-infrared laser therapy**).

A special feature of the RIKTA® – ESMIL® (2A) device is the possibility of **TENS\*-therapy** (the impact of low-frequency pulsed electric current). As a rule, TENS therapy is effectively used for relief of the acute and chronic pain.

 TENS-therapy can be carried out both in conjunction with photo-magnetic-infrared laser therapy, and separately from it.

\* TENS – Transcutaneous Electrical Nerve Stimulation.

## GENERAL INFORMATION

- The present operating manual (further – the Manual) is the document that certifies basic parameters and technical specifications of the electro-photo-magnetic-infrared laser therapy device RIKTA® – ESMIL® (2A), guaranteed by the Manufacturer, and is intended for acquaintance with the operating rules of the device. The manual is also a guide for technical service, transportation and storage.

**The Operating manual should be kept during all the Device service life.**

- It is necessary to read the Manual carefully before using the device. Particular attention should be paid to compliance with the safety precautions listed on page 28, as well as to the list of medical contraindications to laser therapy or TENS-therapy.
- Methods of carrying out treatment procedures are stated in the “Methodical manual on application of the laser therapy device RIKTA® – ESMIL® (2A)” (further, the Methodical manual, which is a part of the delivery set of each device.
- The manufacturer reserves the right to make adjustments and amendments to this Manual and the technical documentation at any time without notice.
- The manufacturer reserves the right to make changes in the configuration and design of the device, which will not affect its operation and basic technical specifications.

EN

The Device and its components should be disposed in accordance with general applicable standards and local environment protection rules.

## SAFETY PRECAUTIONS

Using the Device is allowed only after acquaintance with this Manual and the Methodical manual.

### DON'T OPEN UP

- Never open the Device and never try to disassemble or modify its internal components. The Device has no components serviceable by the user. In case of malfunction, immediately stop using the Device and resort to the manufacturer or the service center.

### PLUG AND POWER CORD

- Do not plug the battery charger into mains with wet hands to avoid electric shock.
- Do not pull the battery charger cord when disconnecting the power cord, make sure to hold the plug when disconnecting the power cord.

### WHEN MALFUNCTIONS ARE OBSERVED

- In case of wear or damage of the cord or battery charger plug, as well as in case of strange odor and smoke, immediately turn off the Device, disconnect the power cord from mains and resort to the manufacturer or the service center.

### EXTERNAL EFFECTS

- It is prohibited to spray flammable substances in the vicinity of the Device.
- Avoid exposure of the Device to moisture and excessive heat.
- The device complies with ISO 60601-1-2 standard for electromagnetic compatibility. However, as a precaution it is not recommended to use the Device nearby the devices sensitive to electromagnetic radiation (e.g. mobile phones, medical equipment, computers, etc.).

## WHEN USING THE DEVICE

- The Device falls within Class 1M laser equipment.

**When using the Device, it is recommended to avoid eyes exposure to direct or reflected radiation. Do not stare into the beam and do not make direct observations using optical instruments.**

- Medical personnel are recommended to wear goggles while operating the Device (available on request). It is allowed not to use goggles when carrying out a single treatment procedure on condition that the Device aperture is tightly pressed to a patient's body.


If the treatment is carried out in a patient's face area with the help of light guide nozzles, the goggles are also recommended for the patient.

- Start the laser therapy procedure only after putting the Device aperture on the treatment zone.

## DEVICE CARE

- Before and after a treatment procedure, make sure to wipe the Device with a cloth of gauze or coarse calico soaked in a treatment tool disinfectant based on 3% hydrogen peroxide solution. Avoid Device contact with alcohol-containing media!
- The Device disinfection may be carried out only when the Device is turned off and disconnected from mains.

## OPERATION CONDITIONS

The Device is designed for operation at a temperature from +10 to +35 °C, atmospheric pressure from 630 to 800 mm Hg, and relative humidity not more than 80%.

## MEDICAL CONTRAINDICATIONS TO LASER THERAPY

**It is STRICTLY FORBIDDEN to carry out laser therapy at home in presence of the following medical contraindications:**

- individual intolerance of influence factors of the device;
- general serious condition of the patient;
- blood diseases with the syndrome of bleeding;
- thrombocytopenia below 60,000 /  $\mu$ l;

- acute feverish conditions;
- active tuberculosis;
- decompensated condition of cardiovascular, excretory, respiratory and endocrine systems;
- deep vein thrombosis;
- malignant neoplasm;
- benign neoplasms with a tendency to grow and become malignant;
- presence of implanted pacemaker;
- all terms of pregnancy;
- period 2 days before and during menstruation;
- children under 5 years of age.

**While treating the apex cardiac beat zone or other cardiac projection, 5 Hz is the only frequency permitted to be used.**

### MEDICAL CONTRAINDICATIONS TO TENS-THERAPY

**It is STRICTLY FORBIDDEN to carry out TENS therapy in presence of the following medical contraindications:**

- individual intolerance of electric current and metal electrodes;
- presence of implanted pacemaker (artificial pacemaker);
- status epilepticus and convulsive disorder
- tumors of any etiology and localization;
- acute febrile condition;
- deep venous thrombosis;
- condition of acute mental, alcohol or drug abuse;
- diseases and injuries of the skin in areas of electrode application;
- children under 14 years of age.

EN

**It is STRICTLY FORBIDDEN to carry out TENS-therapy on any projection of the heart, as well as on throat, perineum, zones of localization of large arteries.**


**To avoid burns during TENS therapy, increase the voltage on the electrodes carefully, by briefly pressing the «TENS +» button (see page 38)**

## DELIVERY SET

No	Description	Quantity
<b>K1</b>	RIKTA® – ESMIL® (2A) device	1 pc
<b>K2</b>	Protective light-guide nozzle (at the time of delivery, screwed in device aperture)	1 pc
<b>K3</b>	Charger	1 pc
<b>K4</b>	Operating manual	1 pc
<b>K5</b>	Methodical manual for using RIKTA® – ESMIL® (2A) device	1 pc
<b>K6</b>	Device storage bag	1 pc

# DEVICE

## K1 – RIKTA®- ESMIL® (2A) device

- 1 – pulse frequency LED indicators:  
5 Hz, 50 Hz, 1000 Hz, VAR – variable frequency;
- 2 – button for setting the required pulse repetition frequency /  
during device operation – to reduce voltage on the electrodes;
- 3 – symbol of presence of laser radiation;
- 4 – treatment time LED indicators: 1 min, 2 min, 5 min;
- 5 – button for setting treatment time /  
during device operation – to increase voltage on the electrodes;
- 6 – «Battery» – battery charge indicator;
- 7 – «TENS» – TENS indicator;
- 8 – button «START/STOP»;
- 9 – socket for connection of charger;
- 10 – Device aperture;
- 11 – electrodes;
- 12 – nameplate:

- L1 – manufacturer's logo;
- L2 – applied part of BF-type;
- L3 – symbol of the presence of laser radiation;
- L4 – laser safety class;
- L5 – sign of certification of Gosstandart of  
the Russian Federation;
- L6 – manufacturer;
- L7 – output parameters of the power  
adapter for charging battery;
- L8 – polarity for charging battery;
- L9 – serial No.

## K2 – Protective light-guide nozzle

At the time of delivery, a protective light-guide nozzle is screwed in the Device aperture [10]. It performs a dual function: to protect the Device aperture from external influences (e.g., dust and dirt) and to allow conducting therapeutic cosmetic procedures (e.g., smoothing wrinkles).


## K3 – Charger

It is used to charge the built-in battery of the device.

## PREPARATION FOR OPERATION

❗ After staying at low temperature the Device shall be kept at a room temperature in the range of  $+20 \pm 5^\circ\text{C}$  for at least 2 hours before switching it on.

### DEVICE ACTIVATION

If, after switching on the device, only the “5 Hz” and “2 min” indicators light up you can proceed directly to the treatment procedure (see “Conducting the therapy session”, page 36).

If the Device had been off for more than **8 hours** then when the Device is turned on **all frequency indicators [1] and all time indicators [4]** will light up.

At that buttons  and  would not respond to pressing: the Device is standing by waiting on the serviceability check (see “Device Serviceability Check”, page 35).

If the **Battery» indicator [6] lights red** after turning on or while working with the device, it means that the battery is discharged and that the battery must be recharged for normal device operation (see “Battery recharging”, page 40).


## DEVICE SERVICEABILITY CHECK (TESTING)

If the Device had been off for more than 8 hours then when the Device is turned on, the serviceability check standby mode will be activated automatically. At that all frequency [1] indicators and all time [4] indicators will light up.

 and  buttons would not respond to pressing.

### For the Device serviceability check:

 Place the device so that the device aperture [10] is not directed to a light source.

If the indicators «5 Hz» and «2 min» only light up after the beep, **device is functional and ready for use** (see «Conducting the therapy session», page 36).

If all pulse frequency LED indicators [1] and all treatment time LED indicators [4] blink and you hear an intermittent sound signal, it means **the Device is faulty**.

Turn it off immediately and contact the manufacturer or service center for assistance.

## CONDUCTING THE THERAPY SESSION

- i** The sequence of the session, the zones, on which the treatment impact is carried out, frequency of laser radiation and time of exposure are indicated in the Methodical manual or in an individual doctor's prescription.


For each zone, the following types of treatment impacts are available:

- **laser therapy without TENS-therapy** (Frequency is specified in the column "Frequency" of the Methodical Manual, a dash is in the column "TENS")
- **laser therapy in combination with TENS-therapy** (Frequency is specified in the column "Frequency" of the Methodical Manual, TENS intensity is specified in the column "TENS")
- **only TENS-therapy** (a dash is in the column "Frequency" of the Methodical Manual, TENS intensity is specified in the column "TENS").

- i** **It is impossible to carry out TENS-therapy with a protective light-guide nozzle [K2] screwed in the device aperture.**  
Unscrew it by turning it counterclockwise and put it in a safe place.

### For each zone:

- I.** In accordance with the recommendations contained in the Methodical Manual or in an individual doctor's prescription:

- by pressing the button  sequentially, set the required laser frequency: 5 Hz, 50 Hz, 1000 Hz, VAR (variable frequency).

If none of the frequency indicators [1] is illuminated (a dash is in the column «Frequency» of the Methodical Manual), only TENS-therapy is conducted;

- by pressing the button  sequentially, set the exposure time: 1 min, 2 min, 5 min.

II. Set the device aperture [10] to the specified zone of treatment.


III. Press the  button.

Short double beep will sound:

**treatment impact on the zone begins.**

During the treatment procedure, time LED indicators [4] will blink and indicate the time remaining until the end of the exposure:

- more than 2 minutes left – «5 min» indicator blinks;
- 1-2 minutes left – «2 min» indicator blinks;
- less than 1 minute left – «1 min» indicator blinks.

**From the moment of start and till the end of treatment, the buttons**

  **and**  **will work in the mode:**

- «TENS -» (to reduce the voltage of electrodes)
- «TENS +» (to increase the voltage of electrodes).

 **At the beginning of a treatment procedure, the voltage is not supplied on the electrodes.**

- If the Methodical manual indicates that the impact on the zone is carried out without TENS-therapy, after the beep do not press the button  before the end of exposure.
- If the Methodical manual indicates that the impact on the zone is carried out with TENS-therapy, after the beep, press the button  once, to start TENS-therapy. TENS indicator [7] will start to blink. Then adjust the TENS intensity (voltage on electrodes) according to the personal sensations of the patient by short pressing of the  (TENS+) and  (TENS -) buttons.

 **To avoid unpleasant sensations and burns, increase voltage on the electrodes gently by shortly pressing the button  (TENS +).**

*Read more about TENS-therapy on page 38.*

On expiry of the set time, a triple “beep” will sound: **the treatment of the zone is finished.**

** You can stop treatment at any time by pressing  button.**

If you press this button again, the treatment will start again at the previously set frequency and time values.

If it is necessary to carry out TENS-therapy, the voltage level on the electrodes is set again.

## TENS-THERAPY

After starting the treatment and before its completion, you can start TENS-therapy at any time with a short press of the button .

** Only when the electrodes are in contact with patient's skin, the electrodes will be supplied with the voltage.**

Then, briefly press the buttons  (TENS +) and  (TENS -) to adjust the TENS intensity (electrodes voltage).

TENS intensity is adjusted by personal sensations of a patient (see the Methodical manual):

- **minimum TENS intensity (Min.):** the patient feels a subtle tingling of the skin under the electrodes;
- **optimal TENS intensity (Optim.):** the patient feels a nonirritating tingling of the skin under the electrodes, when scanning, a specific sound is heard under the electrodes;
- **maximum TENS intensity (Max.):** the patient feels a distinct tingling and vibration, a specific sound is heard under electrodes, there is resistance to movement of electrodes on skin.

TENS indicator [7] blinks during the TENS-therapy.

The blinking frequency of the TENS indicator depends on electrodes voltage: the higher the voltage level is, the higher the blinking frequency is.

The color of TENS indicator depends on the contact with the patient's skin:

<i>Green</i>	Contact with patient's skin is good.
<i>Blue</i>	Contact with patient's skin is good. The so called "saturation" of the patient's skin occurs under the electrodes (a beep sounds at the moment of "saturation"). If there are no specific instructions from the physician, continue the treatment.
<i>Green/Blue</i>	No contact with the patient's skin.

During the device transfer from one zone to another before the end of procedure, the voltage level on the electrodes is maintained at the set value.

After the time of procedure expires (or, if the procedure was interrupted by pressing the button ) , the impact of the TENS stop.

When you restart the treatment procedure, the TENS intensity should be set again.

## BUILT-IN BATTERY MONITORING

- i** After the Device is switched on and during Device operation, battery charge is checked automatically.


The battery charge indicator “Battery” [6] shows battery charging status.

<i>The indicator does not light</i>	The battery voltage is normal
<i>The indicator blinks red</i>	Battery will last for approximately 10 minutes only. It should be recharged.

In case of low battery the Device will be turned off.

## BATTERY RECHARGING

Connect the charger [K3] to the socket [9]

The battery charge indicator [6] lights red during charging process.

**Operating the device while charging the battery is not possible!**

After the charging process is complete, the indicator [6] lights up green. Unplug the charger from the wall outlet and charger cable from the socket of the Device [9].

## DEVICE TURNING OFF

After you finish using the device, turn it off: press and hold the button  for 3–5 seconds.

All indicators will go out, device will switch off.

-  If the device in the switched-on state is not used for 1 minute (no pressing of buttons, no treatment procedure is performed), the device will turn off automatically.

## TECHNICAL SPECIFICATIONS\*

### Photo-magnetic-infrared laser therapeutic mode:

Radiation wavelength,  $\mu\text{m}$ :

- laser infrared radiation ..... 0,80...0,91
- broadband infrared radiation..... 0,86...0,96
- visible red radiation ..... 0,60...0,67

Laser and broadband infrared radiation pulse repetition frequency, Hz:

- fixed ..... 5, 50, 1000
- variable in the range..... 250...1

Red light pulse repetition frequency, Hz..... 2

Magnetic field induction, mT .....  $35 \pm 10$

Laser radiation pulse power, W .....  $16^{+2}_{-4}$

Broadband infrared radiation average power, mW .....  $60 \pm 30$

Red radiation average power, mW ..... 7

Device aperture area,  $\text{cm}^2$  .....  $4 \pm 0,4$

Exposure time, min. .... 1, 2, 5

### TENS mode:

Amplitude of output voltage pulses on metal electrodes, V..... 0...400

Impulse duration at 0.5 of maximum amplitude,  $\mu\text{s}$ , at least..... 20

Impulse frequency, Hz ..... 60...120

### Other parameters:

The Device power supply from built-in battery

- voltage, V ..... 3,7
- capacity,  $\text{mA} \cdot \text{hour}$  ..... 1800

Charger output voltage, V ..... 5

Charger input voltage.....  $110...240 \text{ V}/50 \text{ Hz}$

Electrical safety ..... protection class II

Applied part ..... type BF

Laser safety by IEC 60825-1-2014 .....class 1M

The device complies with IEC 60601-1-2:2014 standard on electromagnetic compatibility.

Net weight (w/o charger), g, no more, .....250

Service life, years ..... 5

EN

\* For reference

## MAINTENANCE

- To ensure reliable operation of the device, perform its maintenance in a timely manner, guided by the table 1.

Table 1

Maintenance interval	Job content	Requirements
Before each session	Visual inspection without using special tools	No mechanical damage of the device case and the charger
	Check the following:	
	1. Battery charging status according to the battery charge indicator [6]	The battery charge indicator [6] should not light
	2. Presence of radiation at the device aperture [10]	After setting of device aperture [10] on the specified exposure zone and pressing of the button  , red light should be visible.
3. Disinfection (if necessary)	In accordance with the recommendations of the section «SAFETY PRECAUTIONS», page 28	

- During maintenance, please observe the safety precautions listed under “SAFETY PRECAUTIONS”, page 28.
- In case of detection of noncompliance of the device with the requirements specified in table 1, further operation of device is not permitted, the device is to be repaired.

## POSSIBLE MALFUNCTIONS AND TROUBLESHOOTING

The list of possible malfunctions, causes and troubleshooting are listed in table 2.

Table 2

Malfunctions	Causes	Troubleshooting
After pressing the button  none of the indicators [1] and [4] do not light up	The battery is completely discharged or the Device is defective	Charge the battery (see page 40). If the indicators [1], [4] do not light up after charging the battery and turning on the device, repair is carried out by the Manufacturer
Any of the indicator [1] and [4] does not light up when the device is switched on	Trouble of normal operation of the device	Repairs are carried out by the Manufacturer
During the functional test (see page 35) all indicators blink [1], [4], an intermittent beep is heard	One or more types of radiation are absent	Repairs are carried out by the Manufacturer

## PACKAGE, STORAGE AND TRANSPORTATION

- Each Device is supplied with protective light-guide nozzle, which is screwed in Device aperture [**K2**]. Each Device, together with accessories, is packed in a storage bag.
- The Device should be stored in the Manufacturer's package in the heated room at temperature in the range of +10...+35 °C and relative humidity of no more than 80% (at +25 °C). The storage room should be free of dust, mercury vapor, acids and alkalis provoking corrosion.
- The Device should be transported by all types of the closed transport in Manufacturer's package. It is necessary to ensure a stable position of the package during transportation.

## MANUFACTURER'S WARRANTY

- The Manufacturer guarantees that the Device parameters correspond to the requirements of the Technical Specifications TU 9444-007-18537172-2003, within 24 months from the date of sale of the Device to the consumer (if the consumer complies with the conditions of use, transportation and storage). Within the aforesaid period, the manufacturer provides free repairs of the Device.
- The Manufacturer is not responsible for and does not make free repairs of malfunctions caused by a fault of the consumer (for example, due to careless handling, transportation, non-compliance with the rules of use, as well as in case of breaking of the seal)
- This warranty does not apply to any of the following:
  - absence of the present Manual;
  - mechanical damage including the damage during transportation;
  - breaking of the seal.
- The warranty also becomes void if during the warranty period a repair of the faulty Device was carried out by unauthorized persons and organizations.
- In the absence of the sales date and the seller's stamp in the Warranty card, the warranty period shall be calculated from the date of Manufacturing of the Device.
- All issues related to repairs and maintenance should be directed to the Manufacturer or to the technical centers listed on the website [www.rikta.ru](http://www.rikta.ru) or in a separate warranty booklet (if available).

After the end of service life, the Manufacturer is not responsible for the safety and compliance of the product with the functional characteristics stated in this Manual.

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